## ZAA 2300 Automatic Film Applicator

- ▶ Automatic laboratory equipment for accurate and reproducible application of coating materials, adhesives and similar products with almost all film applicators
- ▶ Multifunctional use with reversible, double sided glass plate (glass / printing blanket), easy to turn over or to change
- ▶ Adjustable application area with moveable start and stop positions
- ▶ Reproducible, operator independent results
- ▶ Optional vacuumplate
- ▶ Optional versions equipped with isolation plate for use with heatable vacuum plates or heating plates
- ▶ Easy to handle
- ▶ Reliable results
- ▶ Functional design

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Manual applications may have irregularities caused by uneven pressure or non-uniform drawing speed. The ZAA 2300 Automatic Film Applicator ensures an accurate and reproducible coating.

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### Features

- Double sided glass plate for dual-purpose operation: one side with printing blanket for wire-bar applicators and profile rods, on the other side glass surface for other applications
- Option: vacuumplate for substrates; the vacuum created fixes the substrate evenly in place

### Standard delivery

- 1 automatic film applicator
- 1 spirit level
- 1 reversible, double sided glass plate / printing blanket
- 1 power supply cable
- 1 certificate of manufacturer
- 1 polyester mat

### Options

- ACC121/ACC211 precision-vacuumplate with a series of holes of 1 mm (0.04") and 0.5 mm (0.02") respectively
- ACC039/ACC122 (230 V / 115 V) vacuum pump with hose
- ACC1477.230V/ACC1477.115V (230 V / 115 V) heating plate set
- Precision-vacuumplate with a series of holes of 1 mm (0.04") or 0.5 mm (0.02") (230 V or 115 V), heatable up to 150°C (302°F), those options can only be used together with an automatic film applicator equipped with an isolation plate (2300.H/.FH/.FFH)
- ZTC 2200 Test Charts
- Film applicators with maximum outer width of 300 mm (11.81")

### Handling

- Choose the applicator(s).
- Choose the suitable base (glass, vacuum- or heating plate for gap applicators, printing blanket for wire-bar applicators or profile rods).
- Put the wire-bar applicator or profile rod into the ZAA 2300 Automatic Film Applicator and fix with the weight.
- Position other applicators in front of the weight.
- Adjust the desired application length.
- Fix the desired substrate, e.g. test chart, with the start position clamp.
- Set the desired drawing speed.
- Pour the product to be tested in front of the applicator.
- Start the ZAA 2300 Automatic Film Applicator.
- After the application remove the applicator, the weight and the substrate.

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