## ZGR 2020 - 2024 Grindometers

- ▶ Precision grindometers according to Hegman for determination of fineness of grind of coating materials, printing inks, pastes and similar products
- ▶ Easy to handle
- ▶ Easy to clean

## ZGR 2020.A - 2023.A Grindometers

- ▶ Lightweight versions made of aluminium
- ▶ Particularly suitable for determination of fineness of grind of transparent coating materials, printing inks, pastes and similar products
- ▶ Easy to handle
- ▶ Easy to clean

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The fineness of grind is particularly important during the manufacture of pigmented products. By its determination the efficiency of dispersion processes can be followed and the performance of dispersion devices and methods can be determined.

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**Application areas**

- For very different industries, e.g. for the paint, varnish, printing ink, food, cosmetics, paper and detergents industries
- Laboratory test instrument for quality control and experimental purposes
- Applicable to practically all pigmented systems

**Standard delivery**

- 1 grindometer
- 1 certificate of manufacturer
- 1 scraper
- 1 carrying case

**Options**

- ZRT 2025 Groove Depth Meter
- ACC591 calibration and certification (incl. certificate)

**Handling**

- Pour the product under test at the deepest point into the grooves of the grindometer.
- Place the scraper edgewise and parallel to the width of the grindometer on the deepest point of the grooves and draw it within a few seconds (see the standards) to a point beyond the zero depth of the grooves.
- Within seconds (see the standards) after the completion of the drawdown of the sample, view the grindometer at a right angle to the length of the grooves and at an angle of 20° to 30° to the surface of the grindometer in a light that will make the pattern of the sample in the grooves visible.
- Observe the point where in the grooves particles in a greater number (in accordance with ISO/EN/DIN 5 to 10 particles in an imaginary 3 mm [0.12"] wide band) are first visible and read the value on the scale belonging to that point.
- Afterwards clean the grindometer and scraper.
- Check grindometer periodically by use of a groove depth meter.

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Zehntner GmbH Testing Instruments, Gewerbestrasse 4, CH-4450 Sissach, Switzerland/Schweiz