ZTJ 3020 Temperature Control Jacket

- ▶ Three-legged vessel with spirit level for maintaining consistent temperature during viscosity determination and for holding the cup in a level position
- ▶ Viscosity determination at temperatures different from ambient temperature
- ▶ Applicable to all flow cups corresponding to the standards
- ▶ Easy to handle

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Consistent temperature during viscosity determination ensures meaningfully results.

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Application areas

- For paint and varnish manufacturers, their customers and test laboratories
- For paints, varnishes and similar liquids
- The test can be carried out in the laboratory and during production
- Applicable to all common flow cups corresponding to the standards

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Standard delivery

- 1 temperature control jacket for flow cups
- 1 closed double featured vessel for indirect tempering
- 2 tube accesses to temper liquids
- 1 spirit level for levelling
- 1 certificate of manufacturer

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Options

- ZFC 3010 - 3011 as well as ZFC 3013 Flow Cups
- ACC050 insertion thermometer up to 600°C (1112°F)

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Handling

- Thanks to the adjustable feet and spirit level, the temperature control jacket is easy to level.
- Put the chosen flow cup into the temperature control jacket.
- For tempering you need additionally a circulation thermostat.
- Connect the circulation thermostat by use of two heating / cooling tubes on the tube access of the temperature control jacket so that you have a closed circuit.
- Adjust the required temperature at the circulation thermostat and begin with tempering.
- Start the viscosity determination at the required temperature in accordance to the separate instructions for the flow cup and the corresponding standard.
- After finishing the viscosity determination, firstly remove one heating / cooling tube, cover the tube access with your thumb, remove the second heating / cooling tube and pour the temperature control liquid back into the circulation thermostat.

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